

Protocol for Routine Surgery Sanitation

1. Clean all surfaces with solution of _____ (detergent) _____ & hot water
 - a. Include walls, floors, and ALL kennel surfaces including cage bars & ceiling
 - b. thorough cleaning is required to remove surface dirt and to prevent inactivation of the disinfectant
2. Rinse thoroughly with clean water
3. Spray all surfaces with _____ (disinfectant) _____, to continue cleaning procedure.
4. Wipe down with paper towels
5. Re-spray all surfaces with disinfectant at same concentration
6. Let stand appropriate contact time of _____ (contact time) _____ recommended by manufacturer. (if using Peroxigard, rinsing NOT needed)

7. Rinse with clean water

READY FOR USE